

Pelvic Muscle Trainer

XFT-0010EA



Product Components

The wireless pelvic floor muscle stimulator XFT-0010EA is composed of a main unit, electrode probe, charging cable, and app software.



Download on the
App Store

iOS



Android

The XFT-0010EA wireless pelvic floor muscle stimulator uses the host and app to control the output of electric stimulation from the electrode probe, thereby increasing female pelvic floor muscle strength and relieving urinary incontinence, improving sexual quality of life. It is applicable to the treatment of pelvic floor muscle functional disorders (stress, mixed, and moderate urgency urinary incontinence).

Range of Applications

This product is suitable for female pelvic floor muscle training.

Product Features



Self-contained Design



Wireless Probe



Waterproof & Dustproof



APP Controlled



Electrical Stimulation




Bluetooth Connection

Operation Method

The patient lies down or sits semi-upright, placing the electrode probe (with water-soluble lubricant if needed) into the vagina for full contact with pelvic muscles.



As needed, apply a water-soluble lubricant to the electrode probe surface.



Place the electrode probe at the appropriate depth.